

## Anti-SOX10 antibody [R05-9Y-3] (STJA0033748)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P/ICC/IF    |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                      |
|---------------------|------------------------------------------------------|
| Clonality           | Monoclonal                                           |
| Conjugation         | Unconjugated                                         |
| Purification        | Protein A                                            |
| Dilution Range      | WB 1:500-1:2000<br>IHC-P 1:200-1:1000<br>ICC:IF 1:50 |
| Formulation         | Liquid in PBS, Glycerol and BSA                      |
| Isotype             | IgG                                                  |
| Storage Instruction | Store at -20°C. Avoid freeze/thaw cycles.            |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">6663</a>        |
| Gene Symbol | <a href="#">SOX10</a>       |
| UniProt ID  | <a href="#">SOX10_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|