

**Anti-CD169 antibody [R07-8P-1] (STJA0033754)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-P/ICC/IF/FC
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC-P 1:200-1:1000 ICC:IF 1:50 FC 1:50-1:100
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgG
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">6614</a>
Gene Symbol	<a href="#">SIGLEC1</a>
UniProt ID	<a href="#">SN_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---