

Anti-BAP31 antibody [R01-9V-3] (STJA0033766)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS, Glycerol and BSA
Isotype	IgGk
Storage Instruction	Store at -20°C. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	10134
Gene Symbol	BCAP31
UniProt ID	BAP31_HUMAN

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



BAP31 antibody western blot HeLa 293 4T1 rat liver STJA0033766