

## Anti-SOX9 antibody [R06-7F6] (STJA0033897)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P/IP        |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES



|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.3 mg/mL                                                                                |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IP 1:20                                             |
| Formulation         | 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA  |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">6662</a>       |
| Gene Symbol | <a href="#">SOX9</a>       |
| UniProt ID  | <a href="#">SOX9_HUMAN</a> |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SOX9 antibody western blot 3T3 HeLa SOX9 antibody immunohistochemistry human  
 STJA0033897 breast cancer STJA0033897