

## Anti-beta III Tubulin antibody [R02-4G7] (STJA0034079)

### GENERAL INFORMATION

|               |                          |
|---------------|--------------------------|
| Product Type  | Primary antibodies       |
| Applications  | WB/IHC-F/IHC-P/ICC/IF/IP |
| Host / Source | Rabbit                   |
| Reactivity    | Human/Mouse/Rat          |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.3 mg/mL                                                                                |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IF 1:50-1:200<br>IP 1:20                            |
| Formulation         | 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA  |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">10381</a>      |
| Gene Symbol | <a href="#">TUBB3</a>      |
| UniProt ID  | <a href="#">TBB3_HUMAN</a> |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




beta III Tubulin antibody western blot mouse brain  
STJA0034079



beta III Tubulin antibody western blot rat brain  
HeLa STJA0034079



beta III tubulin antibody immunohistochemistry  
human cholangiocarcinoma STJA0034079



beta III Tubulin antibody immunocytochemistry  
HeLa STJA0034079