

Anti-Cadherin antibody [R04-5A8] (STJA0034688)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Hamster

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID
[999](#)
[1001](#)
[1000](#)
[CDH1](#)
[CDH3](#)
[CADH1_HUMAN](#)
[CADH3_HUMAN](#)
[CADH2_HUMAN](#)

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




pan Cadherin antibody western blot mouse heart brain STJA0034688



pan Cadherin antibody western blot HeLa CHO-K1 C6 rat brain STJA0034688



pan-Cadherin antibody immunohistochemistry human cholangiocarcinoma STJA0034688



pan Cadherin antibody immunocytochemistry MCF-7 STJA0034688