

Anti-Phospho-PKA RII alpha-Ser99 antibody [R05-7G2] (STJA0034728)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IP 1:20
Formulation	50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5576
Gene Symbol	PRKAR2A
UniProt ID	KAP2_HUMAN

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Phospho-PKA RII alpha antibody western blot 3T3 Raw2647 STJA0034728



Phospho-PKA RIIalpha Ser99 antibody western blot rat brain C6 STJA0034728



PKA R2 Phospho-Ser99 antibody immunohistochemistry human tonsil STJA0034728