

Anti-RPS6 antibody [R05-3A4] (STJA0034824)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	2.2 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	6194
Gene Symbol	RPS6
UniProt ID	RS6_HUMAN

ADDITIONAL INFORMATION



Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




RPS6 antibody immunohistochemistry human breast cancer STJA0034824



S6 Ribosomal Protein antibody western blot Ramos Jurkat C6 3T3 HeLa STJA0034824



S6 Ribosomal Protein antibody immunoprecipitation HeLa cells STJA0034824



S6 Ribosomal Protein antibody western blot CHO-K1 Raji Raw2647 COS7 STJA0034824