

Anti-RUNX antibody [R07-9C7] (STJA0034832)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	860 861 864 RUNX2 RUNX1 RUNX2_HUMAN RUNX1_HUMAN RUNX3_HUMAN
---------	--

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




RUNX antibody immunocytochemistry Jurkat cells STJA0034832



RUNX antibody immunohistochemistry mouse colon paraffin-embedded STJA0034832



RUNX1/2/3 antibody western blot RAW2647 lysates STJA0034832



RUNX1/2/3 antibody western blot Jurkat rat brain C6 CHO-K1 STJA0034832