

Anti-Phospho-c Myb-Ser11 antibody [R04-2D3] (STJA0034999)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IP/ChIP
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IP 1:20 ChIP 1:20
Formulation	50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	4602
Gene Symbol	MYB
UniProt ID	MYB_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-c-Myb Ser11 antibody western blot
 Jurkat STJA0034999



vMyb/c-Myb Phospho-Ser11 antibody
 immunohistochemistry human tonsil
 STJA0034999