

## Anti-MST1 antibody [R04-1D6] (STJA0035025)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.3 mg/mL                                                                                |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100                                                        |
| Formulation         | 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA  |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">4485</a>       |
| Gene Symbol | <a href="#">MST1</a>       |
| UniProt ID  | <a href="#">HGFL_HUMAN</a> |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



MST1 antibody western blot C2C12 lysates  
STJA0035025



MST1 antibody western blot HeLa C6 rat brain  
STJA0035025



MST1 antibody immunohistochemistry human colon cancer STJA0035025