

Anti-Phospho-delta 1 Catenin-Thr916 antibody [R07-5A6] (STJA0035059)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	1500
Gene Symbol	CTNND1
UniProt ID	CTND1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-delta 1 Catenin antibody western blot
MCF-7 U251 C2C12 STJA0035059



delta 1 Catenin/CAS Phospho-Thr916 antibody
immunocytochemistry HeLa STJA0035059



Phospho-delta 1 Catenin/CAS antibody
immunohistochemistry human lung cancer
STJA0035059