

## Anti-Phospho-eIF4B-Ser406 antibody [R05-4B2] (STJA0035063)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.3 mg/mL                                                                                |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100                                                        |
| Formulation         | 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA  |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">1975</a>       |
| Gene Symbol | <a href="#">EIF4B</a>      |
| UniProt ID  | <a href="#">IF4B_HUMAN</a> |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-eIF4B Ser406 antibody western blot mouse heart STJA0035063



eIF4B Phospho-Ser406 antibody western blot A549 HL-60 U2OS C6 STJA0035063



eIF4B Phospho-Ser406 antibody immunohistochemistry human lung cancer STJA0035063