

Anti-Elongation factor 2 antibody [R07-2C0] (STJA0035189)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	1938
Gene Symbol	EEF2
UniProt ID	EF2_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Elongation factor 2 antibody EEF2 antibody immunocytochemistry HeLa STJA0035189



EEF2 antibody western blot 3T3 MEF EEF2 antibody STJA0035189



EEF2 antibody western blot A549 U2OS C6 EEF2 antibody STJA0035189



EEF2 antibody immunohistochemistry human breast cancer STJA0035189