

## Anti-RanGAP1 antibody [R03-1F8] (STJA0035347)

### GENERAL INFORMATION

|               |                          |
|---------------|--------------------------|
| Product Type  | Primary antibodies       |
| Applications  | WB/IHC-F/IHC-P/ICC/IF/IP |
| Host / Source | Rabbit                   |
| Reactivity    | Human/Mouse/Rat          |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Concentration       | 0.3 mg/mL                                                                                |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Purified                                                                        |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IF 1:50-1:200<br>IP 1:20                            |
| Formulation         | 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA  |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION


|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">5905</a>        |
| Gene Symbol | <a href="#">RANGAP1</a>     |
| UniProt ID  | <a href="#">RAGP1_HUMAN</a> |

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




RanGAP1 antibody western blot 3T3 cells STJA0035347



RanGAP1 antibody western blot HeLa A549 HL-60 U2OS C6 STJA0035347



RanGAP1 antibody immunocytochemistry HeLa cells STJA0035347



RanGAP1 antibody immunohistochemistry human tonsil STJA0035347