

Anti-Ubiquitin Conjugating Enzyme E2 C antibody [R08-4F5] (STJA0035393)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-F/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.3 mg/mL
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	11065
Gene Symbol	UBE2C
UniProt ID	UBE2C_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




UBE2C antibody immunocytochemistry HeLa cells STJA0035393



UBE2C antibody western blot PC-12 Raw2647 3T3 STJA0035393



UBE2C antibody western blot A549 HL-60 U2OS STJA0035393



UBE2C antibody immunohistochemistry human tonsil STJA0035393