

## Anti-V5-Tag antibody [R09-5J4] (STJA0035426)

### GENERAL INFORMATION


|               |                     |
|---------------|---------------------|
| Product Type  | Primary antibodies  |
| Applications  | WB                  |
| Host / Source | Rabbit              |
| Reactivity    | Species-independent |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | Unconjugated   |
| Purification        | Affinity Purified  |
| Dilution Range      | WB 1:1000-1:20000  |
| Formulation         | 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA  |
| Isotype             | IgG  |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



V5 tag antibody western blot fusion protein  
 STJA0035426