

## Anti-Acetyl-Histone H3 antibody [R49-6E-4] (STJA0035496)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:100-1:200 IF 1:50-1:200
Formulation	Liquid in PBS, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">440093</a> <a href="#">3020/3021</a>
Gene Symbol	<a href="#">H3-5</a> <a href="#">H3-3A.H3-3B</a> <a href="#">H3C_HUMAN</a> <a href="#">H33_HUMAN</a> <a href="#">H32_HUMAN</a> <

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Acetyl-Histone H3 Lys27 antibody western blot  
HeLa cells STJA0035496