

## Anti-NDUFAB1 antibody [R34-7I-8] (STJA0035603)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P/IP/FC     |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Chromatography                                                                  |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:100-1:200<br>IP 1:50<br>FC 1:50                                 |
| Formulation         | Liquid in PBS, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                  |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">4706</a>       |
| Gene Symbol | <a href="#">NDUFAB1</a>    |
| UniProt ID  | <a href="#">ACPM_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



NDUFAB1 antibody western blot A431 cells  
STJA0035603