

## Anti-Ceruloplasmin antibody [R18-3Y-7] (STJA0035646)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Rabbit             |
| Reactivity    | Human/Rat          |

### PRODUCT PROPERTIES


|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                               |
| Conjugation         | Unconjugated                                                                             |
| Purification        | Affinity Chromatography                                                                  |
| Dilution Range      | WB 1:500-1:1000                                                                          |
| Formulation         | Liquid in PBS, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                  |
| Isotype             | IgG                                                                                      |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">1356</a>       |
| Gene Symbol | <a href="#">CP</a>         |
| UniProt ID  | <a href="#">CERU_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Ceruloplasmin antibody western blot rat liver  
 STJA0035646