

## Anti-Phospho-CaMKII-Thr286 antibody [R95-4J-2] (STJA0035695)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P
Host / Source	Rabbit
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:100-1:200
Formulation	Liquid in PBS, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

### TARGET INFORMATION

Gene ID  
[25400](#)  
[815](#)  
[12322](#)  
[Camk2a](#)  
[CAMK2A](#)  
[KCC2A\\_RAT](#)  
[KCC2A\\_HUMAN](#)  
[KCC2A\\_MOUSE](#)

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-CaMKII alpha antibody western blot mouse brain STJA0035695