

Anti-PKM2 antibody (STJA0035704)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IHC-IF/FC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB-1:1000-1:2000 IHC-P-1:50-1:200 ICC/IF-1:50-1:200 FC-1:50
Formulation	Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	5315
Gene Symbol	PKM
UniProt ID	KPYM_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



PKM2 antibody western blot COS7 C6 3T3 HeLa
STJA0035704