

Anti-beta Actin antibody (STJA0035722)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB-1:5000-1:10000 IHC-1:50-1:200 IF-1:100-1:200
Formulation	Liquid in 1xPBS (pH 7.4) , 150mM NaCl, 50% Glycerol, 0.02% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	60
Gene Symbol	ACTB
UniProt ID	ACTB_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




beta Actin antibody western blot 3T3 HeLa C6 STJA0035722



beta actin antibody western blot Jurkat lysates STJA0035722



beta Actin antibody immunohistochemistry human heart mouse testis bladder cancer STJA0035722



beta Actin antibody immunofluorescence HeLa cells STJA0035722