

Anti-Estrogen Receptor alpha antibody [R09-7G3] (STJA0035748)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/ChIP/FC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Dog

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 FC 1:50-1:100 ChIP 1:20
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	2099
Gene Symbol	ESR1
UniProt ID	ESR1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ER alpha antibody western blot MCF-7 T47D
STJA0035748



ER alpha antibody immunohistochemistry human
cervix carcinoma STJA0035748



Estrogen Receptor alpha antibody
immunofluorescence MCF-7 STJA0035748