

Anti-Phospho-PKR-Thr446 antibody [R03-5B4] (STJA0035856)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IP/ChIP
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IP 1:50 ChIP 1:20
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	5610
Gene Symbol	EIF2AK2
UniProt ID	E2AK2_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Phospho-PKR antibody western blot HeLa lysates
STJA0035856



Phospho-PKR immunohistochemistry
STJA0035856

T446 human antibody colon