

Anti-Cyclin D1 antibody [R09-7J5] (STJA0035915)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/IP
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:50
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	595
Gene Symbol	CCND1
UniProt ID	CCND1_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Cyclin D1 antibody western blot NIH 3T3 C6 lysates STJA0035915



Cyclin D1 antibody western blot MCF-7 LNCaP lysates STJA0035915



Cyclin D1 antibody immunohistochemistry human bladder STJA0035915



Cyclin D1 antibody immunofluorescence MCF-7 cells STJA0035915