

## Anti-CREB antibody [R04-4A8] (STJA0035922)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IHC-P/ICC/IF/IP/FC |
| Host / Source | Rabbit                |
| Reactivity    | Human/Mouse/Rat       |

### PRODUCT PROPERTIES


|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                        |
| Conjugation         | Unconjugated                                                                                      |
| Purification        | Affinity Chromatography                                                                           |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IF 1:50-1:200<br>IP 1:50<br>FC 1:50-1:100                    |
| Formulation         | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Isotype             | IgG                                                                                               |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |

### TARGET INFORMATION


|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">1385</a>        |
| Gene Symbol | <a href="#">CREB1</a>       |
| UniProt ID  | <a href="#">CREB1_HUMAN</a> |

### ADDITIONAL INFORMATION


**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CREB antibody western blot NIH/3T3 lysates  
STJA0035922



CREB antibody immunohistochemistry mouse  
brain STJA0035922



CREB antibody immunofluorescence HeLa cells  
STJA0035922