

Anti-Phospho-JNK-Thr183 antibody [R04-9D2] (STJA0035957)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/IF/IP/FC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:50 FC 1:50-1:100
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID
[5601](#)
[5602](#)
[5599](#)
[MAPK9](#)
[MAPK10](#)
[MK09_HUMAN](#)
[MK10_HUMAN](#)
[MK08_HUMAN](#)

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Phospho-JNK antibody western blot
 HeLa C6 UV STJA0035957



Phospho-JNK antibody western blot
 NIH/3T3 STJA0035957



Phospho-JNK1/2/3 immunohistochemistry
 human breast tissue STJA0035957



Phospho-JNK antibody immunofluorescence
 NIH/3T3 Anisomycin STJA0035957