

Anti-Tyrosine Hydroxylase antibody [R06-1A8] (STJA0035961)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ICC/IF/FC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 FC 1:50-1:100
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	7054
Gene Symbol	TH
UniProt ID	TY3H_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Tyrosine 3 Monooxygenase antibody western blot rat mouse brain STJA0035961



Tyrosine Hydroxylase antibody western blot PC-3 cells STJA0035961



Tyrosine Hydroxylase antibody immunohistochemistry human kidney STJA0035961



Tyrosine Hydroxylase antibody immunofluorescence PC-12 cells STJA0035961