

Anti-CD3D antibody [R03-2H9] (STJA0036216)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IHC-IF/IP/FC
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200 IP 1:50 FC 1:50-1:100
Formulation	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION


Gene ID	915
Gene Symbol	CD3D
UniProt ID	CD3D_HUMAN

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CD3D antibody western blot Jurkat lysates
STJA0036216



CD3D antibody immunohistochemistry human spleen
STJA0036216



CD3D antibody immunofluorescence Jurkat cells
STJA0036216