

## Anti-VPS24 antibody [R04-9R6] (STJA0037220)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC-P/IP/FC     |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                        |
| Conjugation         | Unconjugated                                                                                      |
| Purification        | Affinity Purified                                                                                 |
| Dilution Range      | WB 1:500-1:1000<br>IHC 1:50-1:100<br>IP 1:50<br>FC 1:50-1:100                                     |
| Formulation         | Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Isotype             | IgG                                                                                               |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |

### TARGET INFORMATION

|             |                                 |
|-------------|---------------------------------|
| Gene ID     | <a href="#">100526767/51652</a> |
| Gene Symbol | <a href="#">CHMP3</a>           |
| UniProt ID  | <a href="#">CHMP3_HUMAN</a>     |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



VPS24 antibody western blot HepG2 lysates  
STJA0037220