

## Anti-Estrogen Related Receptor alpha antibody (STJA0037543)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/FC/IP           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                        |
| Conjugation         | Unconjugated                                                                                      |
| Purification        | Affinity Purified                                                                                 |
| Dilution Range      | WB-1:1000-1:5000<br>FC-1:20-1:50<br>IP-1:20                                                       |
| Formulation         | Liquid in 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA. |
| Isotype             | IgG                                                                                               |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">2101</a>       |
| Gene Symbol | <a href="#">ESRRA</a>      |
| UniProt ID  | <a href="#">ERR1_HUMAN</a> |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ERR alpha antibody western blot HeLa lysates  
 STJA0037543