

## Anti-mCherry-Tag antibody (STJA0037573)

### GENERAL INFORMATION

|               |                       |
|---------------|-----------------------|
| Product Type  | Primary antibodies    |
| Applications  | WB/IHC-P/IHC-IF/FC/IP |
| Host / Source | Rabbit                |
| Reactivity    | Species-independent   |

### PRODUCT PROPERTIES

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                        |
| Conjugation         | Unconjugated                                                                                      |
| Purification        | Affinity Purified                                                                                 |
| Dilution Range      | WB-1:500-1:1000<br>IHC-P-1:50-1:100<br>ICC/IF-1:50-1:100<br>FC-1:200-1:500<br>IP-1:20             |
| Formulation         | Liquid in 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA. |
| Isotype             | IgG                                                                                               |
| Storage Instruction | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|