

Anti-GSDMD antibody (STJA0037577)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/FC/IP
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB-1:500-1:1000 IHC-1:200-1:1000 FC-1:200-1:500 IP-1:20-1:50
Formulation	Liquid in 50mM Tris-Glycine (pH 7.4) , 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	69146
Gene Symbol	Gsdmd
UniProt ID	GSDMD_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---