

Anti-Recombinant-chimeric-CD72-Nb-Fc antibody (STJA0038340)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	293F Cells
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	ELISA 1:5000-100000
Formulation	Phosphate-buffered solution
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	971
Gene Symbol	CD72
UniProt ID	CD72_HUMAN
Specificity	This recombinant monoclonal antibody can detects endogenous levels of CD155 protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---