

Anti-Recombinant-chimeric-CD94-Nb-Fc antibody (STJA0038469)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	293F Cells
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	ELISA 1:5000-100000
Formulation	Phosphate-buffered solution
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	This recombinant monoclonal antibody can detect endogenous levels of CD276/B7H3 protein.
-------------	--

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ELISA: STJA0038469