

## Anti-Recombinant-chimeric-CD166-Nb-Fc antibody (STJA0038597)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	293F Cells
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	ELISA 1:5000-100000 FC 1-2 Mu g/Test
Formulation	Phosphate-buffered solution
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">214</a>
Gene Symbol	<a href="#">ALCAM</a>
UniProt ID	<a href="#">CD166_HUMAN</a>
Specificity	This recombinant monoclonal antibody can detects endogenous levels of CD45 protein.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ELISA: STJA0038597