

Anti-Recombinant-chimeric-CD314-Nb-Fc antibody (STJA0038598)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA |
| Host / Source | 293F Cells |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | Lot specific |
| Conjugation | Unconjugated |
| Purification | Recombinant Expression and Affinity purified |
| Dilution Range | ELISA 1:5000-100000 |
| Formulation | Phosphate-buffered solution |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|--|
| Specificity | This recombinant monoclonal antibody can detect endogenous levels of CD45 protein. |
|-------------|--|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|