

## Anti-Recombinant-chimeric-CD82-Nb-Fc antibody (STJA0038608)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA              |
| Host / Source | 293F Cells         |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                     |
| Concentration       | Lot specific                                                                                   |
| Conjugation         | Unconjugated                                                                                   |
| Purification        | Recombinant Expression and Affinity purified                                                   |
| Dilution Range      | ELISA 1:5000-100000                                                                            |
| Formulation         | Phosphate-buffered solution                                                                    |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|             |                                                                                    |
|-------------|------------------------------------------------------------------------------------|
| Gene ID     | <a href="#">3732</a>                                                               |
| Gene Symbol | <a href="#">CD82</a>                                                               |
| UniProt ID  | <a href="#">CD82_HUMAN</a>                                                         |
| Specificity | This recombinant monoclonal antibody can detect endogenous levels of CD48 protein. |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|