

Anti-Recombinant-chimeric-CD96-Nb-Fc antibody (STJA0038630)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	293F Cells
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	ELISA 1:5000-100000
Formulation	Phosphate-buffered solution
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10225
Gene Symbol	CD96
UniProt ID	TACT_HUMAN
Specificity	This recombinant monoclonal antibody can detect endogenous levels of CD54 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ELISA: STJA0038630