

Anti-Recombinant-chimeric-CD55-Nb-Fc antibody (STJA0038639)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	293F Cells
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	ELISA 1:5000-100000 FC 1-2 Mu g/Test
Formulation	Phosphate-buffered solution
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1604
Gene Symbol	CD55
UniProt ID	DAF_HUMAN
Specificity	This recombinant monoclonal antibody can detects endogenous levels of CD56 protein.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.