

Anti-Recombinant-chimeric-CD16-Nb-Fc antibody (STJA0038648)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/FC
	Host / Source	293F Cells
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Recombinant Expression and Affinity purified
	Dilution Range	ELISA 1:5000-100000 FC 1-2 Mu g/Test
	Formulation	Phosphate-buffered solution
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Specificity	This recombinant monoclonal antibody can detects endogenous levels of CD62L protein.
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.