

Anti-Recombinant-chimeric-CD74-Nb-Fc antibody (STJA0038654)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	293F Cells
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	ELISA 1:5000-100000
Formulation	Phosphate-buffered solution
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	972
Gene Symbol	CD74
UniProt ID	HG2A_HUMAN
Specificity	This recombinant monoclonal antibody can detect endogenous levels of CD62L protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--