

**Anti-Recombinant-chimeric-CD44-Nb-Fc antibody (STJA0038748)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	293F Cells
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	ELISA 1:5000-100000 FC 1-2 Mu g/Test
Formulation	Phosphate-buffered solution
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">960</a>
Gene Symbol	<a href="#">CD44</a>
UniProt ID	<a href="#">CD44_HUMAN</a>
Specificity	This recombinant monoclonal antibody can detects endogenous levels of CD8 protein.

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---