

**Anti-Recombinant-chimeric-CD79b-Nb-Fc antibody (STJA0038848)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA              |
| Host / Source | 293F Cells         |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                               |
|---------------------|-----------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                    |
| Concentration       | Lot specific                                                                                  |
| Conjugation         | Unconjugated                                                                                  |
| Purification        | Recombinant Expression and Affinity purified                                                  |
| Dilution Range      | ELISA 1:5000-100000                                                                           |
| Formulation         | Phosphate-buffered solution                                                                   |
| Storage Instruction | Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

**TARGET INFORMATION**

|             |                                                                                     |
|-------------|-------------------------------------------------------------------------------------|
| Gene ID     | <a href="#">974</a>                                                                 |
| Gene Symbol | <a href="#">CD79B</a>                                                               |
| UniProt ID  | <a href="#">CD79B_HUMAN</a>                                                         |
| Specificity | This recombinant monoclonal antibody can detects endogenous levels of IL-4 protein. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|