

Anti-Recombinant-CD195/CCR5 antibody [PTR1380] (STJA0038947)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	ELISA 1:5000-100000 FC 1-2 Mu g/Test
Formulation	Phosphate-buffered solution
Isotype	IgG1k
Molecular Weight	Observed: 40kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	This recombinant monoclonal antibody can detect endogenous levels of CD195 protein.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--