

Anti-Glycophorin A antibody (1-100) [ABT-GYPA] (STJA0039004)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/IF/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	IHC 1:200-1000 WB 1:500-2000 IF 1:100-500 ELISA 1:1000-5000
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG1k
Molecular Weight	Calculated: 16kD Observed: 37kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Immunogen Region	1-100
Specificity	The antibody can specifically recognize human Glycophorin A protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western Blot (WB): Whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with STJA0039004




Western Blot (WB): Various whole cell lysates were separated by 15% SDS-PAGE, and the membrane was blotted with STJA0039004



Immunohistochemistry (IHC): Human acute myeloid leukemia tissue was stained with Anti-Glycophorin A (ABT-GYPA) antibody. STJA0039004



Immunohistochemistry - Paraffin (IHC-P): Human appendix tissue was stained with anti-Glycophorin A (ABT-GYPA) antibody. STJA0039004



Immunohistochemistry - Paraffin (IHC-P): Human tonsil tissue was stained with Anti-Glycophorin A (ABT-GYPA) Antibody. STJA0039004