

Anti-HDAC6 antibody (1100-1215) [ABT-HDAC6] (STJA0039006)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/IF/ELISA
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	IHC 1:200-1000 WB 1:500-2000 IF 1:100-500 ELISA 1:1000-5000
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG1k
Molecular Weight	Calculated: 131kD Observed: 130180kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10013
Gene Symbol	HDAC6
UniProt ID	HDAC6_HUMAN
Immunogen Region	1100-1215
Specificity	The antibody can specifically recognize human HDAC6 protein.

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western Blot (WB): Whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-HDAC6 (STJA0039006)



Immunohistochemistry - Paraffin (IHC-P): Human testis tissue was stained with Anti-HDAC6 (ABT-HDAC6) Antibody. STJA0039006