

Anti-Nucleophosmin antibody (100-200) [ABT486] (STJA0039082)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/IF/ELISA
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	IHC 1:200-1000 WB 1:500-2000 IF 1:100-500 ELISA 1:1000-5000
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG2ak
Molecular Weight	Calculated: 32kD Observed: 33kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Immunogen Region	100-200
Specificity	The antibody can specifically recognize human Nucleophosmin protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western Blot (WB): A431 whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with STJA0039082



Immunohistochemistry - Paraffin (IHC-P): Human colon carcinoma tissue was stained with Anti-Nucleophosmin (ABT486) STJA0039082



Immunohistochemistry - Paraffin (IHC-P): Human hepatocellular carcinoma tissue was stained with Anti-Nucleophosmin STJA0039082