

Anti-SALL4 antibody (750-850) [PT2242] (STJA0039096)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/WB/IF/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G
Dilution Range	IHC 1:200-1000 WB 1:500-2000 IF 1:100-500 ELISA 1:1000-5000
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgG2bk
Molecular Weight	Calculated: 113kD Observed: 130kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	57167
Gene Symbol	SALL4
UniProt ID	SALL4_HUMAN
Immunogen Region	750-850
Specificity	This antibody detects endogenous levels of SALL4 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western Blot (WB): Various whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with STJA0039096



Western blot analysis of SALL4 Antibody at 1:1000 dilution. STJA0039096



Immunohistochemical analysis of paraffin-embedded human seminoma. 1, SALL4 antibody was diluted at 1:200 (4°C STJA0039096